Caution: For Laboratory Use. A product for research purposes only.

Hydrogen Peroxide Luciferin

Description:

Hydrogen peroxide (H_2O_2) is a component of cell signaling pathways that are necessary for the growth, development, and fitness of living organisms. Hydrogen Peroxide Luciferin is a chemoselective bioluminescent probe for the real-time detection of H_2O_2 in cell culture and living animals. Hydrogen Peroxide Luciferin is a boronic acid-caged firefly luciferin molecule that selectively reacts with H_2O_2 to release firefly luciferin, which triggers a bioluminescent response in the presence of firefly luciferase. The high sensitivity and selectivity of Hydrogen Peroxide Luciferin for H_2O_2 , combined with the favorable properties of bioluminescence for in vivo imaging, afford a unique technology for real-time detection of basal levels of H_2O_2 generated in healthy, living mice.



SwissLumix Sarl EPFL Technology Park, Batiment C Lausanne, Switzerland T: +41 78 857 1660 www.swisslumix.com Selective and concentration-dependent bioluminescent detection of H₂O₂ by hydrogen peroxide luciferin probe (HPL). (A) Total bioluminescent signal, integrated over 45 min, from PCL-1 (5 μ M) alone (light gray bars) or incubated with various ROS (100 μ M) or H₂O₂ (100 μ M) and catalase (0.4 mg/mL) for 5, 20, 40, or 60 min. Signals were normalized to signal from HPL in the absence of any ROS. (B) Total bioluminescent signal, integrated over 45 min, from 5 μ M HPL incubated for 1 h with increasing concentrations of H_2O_2 (0–250 μ M). Signals were normalized to signal from HPL in the absence of H_2O_2 . To quantify free luciferin formation in A and B, 100 µg/mL luciferase in 50 mM Tris buffer with 10 mM MgCl₂, 0.1 mM ZnCl₂, and 2 mM ATP (pH 7.4) was added to the HPL plus ROS solutions.

Bioluminescent signal from HPL in FVB-luc⁺ mice. (A) Representative image (30 min postinjection) for mice injected with HPL (i.p., 0.5 µmol in 50 µL of 1:1 DMSO:PBS) immediately prior to injection of H₂O₂ (i.p., 0, 0.37, 1.5, 6, or 24 mM, left to right, in 100 µL of PBS). (*B*) Total photon flux, integrated over 1 h, for mice injected with HPL \pm H₂O₂. H₂O₂ concentrations represent the H₂O₂ concentration in the i.p. cavity based on a total injection volume of 150 µL. Statistical analyses were performed with a two-tailed Student's *t* test. ***P* < 0.005 (*n* = 5) and error bars are \pm SD. (*C*) Representative image (12 min postinjection) for mice injected with HPL (i.p., 0.5 µmol in 50 µL of 1:1 DMSO:PBS) immediately following NAC (i.p., 0 or 0.2 mg in 100 µL PBS). (*D*) Total photon flux, integrated over 1 h, for mice injected with HPL \pm NAC.



Caution: For Laboratory Use. A product for research purposes only.

Hydrogen Peroxide Luciferin

Imaging and Applications:

- Imaging of hydrogen peroxide fluxes in cell culture and in living animals *in vivo*
- High sensitivity and low background to noise ratio
- Recommended imaging time is 10-60 min post injection of the probe
- Recommended dose is 0.5 uM per mouse injected i.p.

References:

1. Rhee SG. H₂O₂, a necessary evil for cell signaling. Science. 2006;312:1882–1883.

2. Stone JR, Yang S. Hydrogen peroxide: A signaling messenger. Antioxid Redox Signaling. 2006;8:243–270.

3. D'AutrEaux B, Toledano MB. ROS as signaling molecules: Mechanisms that generate specificity in ROS homeostasis. Nat Rev Mol Cell Biol. 2007;8:813–824.

4. Miller EW, Chang CJ. Fluorescent probes for nitric oxide and hydrogen peroxide in cell signaling. Curr Opin Chem Biol. 2007;11:620–625.

5. Winterbourn CC. Reconciling the chemistry and biology of reactive oxygen species. Nat Chem Biol. 2008;4:278–286.

6. Paulsen CE, Carroll KS. Orchestrating redox signaling networks through regulatory cysteine switches. ACS Chem Biol. 2010;5:47–62.

7. Van de Bittnera GC, Dubikovskayaa EA, Bertozzia CR and Changa, CJ. In vivo imaging of hydrogen peroxide production in a murine tumor model with a chemoselective bioluminescent reporter. 21316–21321 | PNAS | December 14, 2010 | vol. 107 | no. 50

SwissLumix Sarl EPFL Technology Park, Batiment C Lausanne, Switzerland T: +41 78 857 1660 www.swisslumix.com

